


Low-cost electrochemical sensing platform for DNA analysis



WELPCR


Supervisor:
Prof. Siddharth Tallur
EE, IIT Bombay



WELSTAT

Student:
Vivekanand Ramrao Dhakane
193079041, *Solid State Devices* (EE-7)
EE, IIT Bombay

WELPCR



Outline

- Rationale
- Objectives
- **WELPCR**
 - Mechanical Assembly
 - Manufacturing of 2D (milling) and 3D (laser cutting) metal parts
 - Electronic Circuit
 - PID design
 - Challenges and Solutions
 - Testing samples using WELPCR
 - Comparison with commercial PCR devices
- **WELSTAT**
 - What is a Potentiostat and Why do we need it?
 - WELSTAT design
 - Selecting AFE for WELSTAT
 - Challenges and Solutions
 - End point PCR using WELPCR+WELSTAT
- Future Improvements

Rationale


The COVID-19 pandemic has highlighted the need of diagnosis especially using the gold standard qRT-PCR
PCR (Polymerase Chain reaction) is a test to tell you if you have an infection

Apart from SARS-CoV-2 (virus that cause COVID-19), PCR can be used to detect other infectious pathogens from environmental samples; blood, saliva and other body fluid samples


1. Sample Collection




2. Sample Processing




3. Amplification & Detection





4. Results



Process-flow for virus detection



Process flow for virus detection using electrochemical sensing



Source: Kumar, M. S., et al. "Electrochemical sensing of SARS-CoV-2 amplicons with PCB electrodes." *Sensors and Actuators B: Chemical* 343 (2021): 130169

Source: Ahuja S, et al. "Longer amplicons provide better sensitivity for electrochemical sensing of viral nucleic acid in water samples using PCB electrodes." *Scientific reports* 12.1 (2022): 1-10.

Virus detection using electrochemical sensor

- The group I am associated with at AIMS lab, EE, IITB have successfully detected SARS-CoV-2 and Phi-6 from water samples using PCB based electrochemical sensing
- While testing the samples, amplification was performed using MiniPCR® thermocycler
- The device costs around \$800-1000+shipping charges (imported from US)
- We faced several maintenance and service issues on the PCR device while performing experiments
- This led to the motivation of designing a PCR device that can be used for virus (or any pathogen) detection using electrochemical sensing



Objectives


- Develop a low-cost, portable, user-friendly PCR device
- Provisions for simultaneous testing of multiple samples (≥ 8 wells)
- Ramp rate ($\geq 1.5^{\circ}\text{C/s}$) similar to commercially available PCR device
- Develop a low-cost Potentiostat to perform electrochemical measurements (Cyclic Voltammetry)




WELPCR




Temperature profile of PCR



WELPCR




WELPCR Block Diagram




Mechanical Assembly



WELPCR: designed using CAD




WELPCR designed using a
Computer aided design tool
Autodesk Fusion 360




WELPCR after
manufacturing

Cross-section view of WELPCR





Manufacturing 3D metal parts



Aluminum block

Milling





Tube holder



Lid bracket



Manufacturing Sheet metal parts



Aluminium sheet after laser cutting

MILESTONE achieved
during MTP Phase-I

*Mechanical Assembly of
WELPCR*




Electronic circuit




Electronic circuit

Circuit assembled on zero PCB to check required components and working of the circuit



H-bridge + micro-controller
+ logic level shifter + display



WELPCR mechanical assembly with sensors

PID Design

(Toughest part of the project)


$$PWM = k_p \times e(t) + k_d \times \frac{de(t)}{dt} + k_i \times \int e(t)dt + \beta$$

β = PWM duty cycle required to maintain constant set-point temperature

- **Trial and error method to tune the PID**
 - ✓ Increase the k_p till we get constant oscillations in steady state
 - ✓ Increase the k_d till the oscillation amplitude is within acceptable range
- **Introducing Integral Controller:** Reduced steady state error but **adds overshoot**
 - ✓ Increase k_i till you get steady state error in acceptable range
- **Reduced k_i to get critically damped system**



PID Tuning



Challenges

1	Poor heat transfer between tube holder and sample
2	Condensation on side walls of tube
3	Noise in temperature sensor
4	Micro-controller out of stock due to chip shortage
5	Accurate temperature sensor calibrations
6	Time required for self-heating of temperature sensor
7	Size of temperature sensor larger than sample size
8	SI2302 MOSFET damage exceeding P_{max} of MOSFET
9	Damage of power supply due to short circuit
10	Spike in output voltage of supply damaging the controller



Challenge: Poor heat transfer


We tried amplifying 503bp fragments of bacteriophage Phi-6 (surrogate for SARS-CoV-2) and failed in getting a band in agarose gel electrophoresis

Reason


Sample inside the test tube could not reach the temperature of tube holder

* Use of 0.5ml tubes

* Tube not covered fully by tube holder



ladder NTC PC



First test using 503bp fragment of Phi6




Tube side walls not fully covered by tube holder

Solutions

Solution 1: Insulation layer


Applied thermal insulation tape covering full body of tube holder



Solution 2: Nichrome heater

Connected nichrome heater to tube holder to heat at higher temperature

Silicone
Nichrome
Ceramic
heat
tape
2mm



Challenge: Condensation

- While performing PCR we observed that no sample was present at the bottom of the tube where actual thermal cycling was being performed

Reason:

- Uncovered part of the tube remains at low temperature causing condensation



Solution


Re-designed tube holder for 0.2ml tubes



Before



After



Challenge: Noise in temperature sensors


- Wrong feedback by sensors to the control system (PID controller)

Reason:

- Peltier is driven by PWM of 12V amplitude, creating Electro Magnetic Interference (EMI) in nearby sensors

Solution:


- Connecting full mechanical assembly to ground of the circuit




Challenges	Solution	Remark
Poor heat transfer between tube holder and sample	Insulation of tube holder	Failed
	Additional nichrome coil heater along with Peltier	Failed
Condensation on side walls of tube	Tube holder designed for smaller tubes covering entire tube sidewalls	Successful (Solved poor heat transfer problem too)
Noise in temperature sensor	Connecting full mechanical assembly to ground of the circuit	Successful
Micro-controller out of stock due to chip shortage	De-soldered microcontroller from old board	Successful
Accurate temperature sensor calibrations	Calculated different equations for each temperature sensor	Successful
Time required for self heating of temperature sensor	Use of smaller 1mm sensor	Successful
Size of temperature sensor larger than sample size	Use of smaller 1mm sensor	Successful
SI2302 MOSFET damage exceeding P_{max} of MOSFET	Use of ON/OFF switch instead of PID controller which operates MOSFET in cutoff or linear region	Successful
Damage of power supply due to short circuit	Use of fuse	Successful
Spike in output voltage of supply damaging the controller	Use of Transient Voltage Suppression(TVS) diode	Successful

Results: Set-up for first successful test

using 503bp fragment of bacteriophage Phi6 (surrogate for SARS-CoV-2)




PCR with live temperature graph on laptop



Mechanical assembly connected to the circuit assembled on zero PCB


ladder NTC PC



Agarose gel electrophoresis result


Results: Second successful test

using 503bp fragment of bacteriophage Phi6 (surrogate for SARS-CoV-2)



Comparison of WELPCR and Bio-Rad

*Comparative image of agarose gel electrophoresis
for 503bp segment of bacteriophage phi6.*




*Estimation of amplified DNA (Nanodrop
Spectrophotometer at 260nm)*


	Thermal cycler	Concentration (ng/μl)
1	Bio-Rad MJ Mini™ Personal thermal Cycler	20.1
2	WELPCR	16.5

Real time sensing using EMStat Potentiostat



Setup for real time testing



Electrode wire connection



Micro-electrode inserted in tube



Real time PCR testing result

Comparison

	WELPCR, EE IITB	OpenPCR (open- source PCR machine design)	miniPCR® mini8	Biorad CFX96 real time PCR	Himedia Wee-32™
Maximum ramp rate	2°C/s	1°C/s	2.4°C/s	2.5°C/s	5°C/s
Reverse transcription (RT-PCR) possible?	Yes	Yes	No	Yes	Yes
Sample capacity	9 wells	16 wells	8 wells	48 wells	32 wells
User interface	Built-in touchscreen and control panel	Software controlled (requires computer)	Mobile app; Software controlled (requires computer)	Built-in control panel and screen	Built-in control panel and screen
Weight	1.9 kg	-	0.45 kg	4.1 kg	3.2 kg
Price	BOM < Rs. 10,000	BOM ~ Rs. 38,500 (USD \$500)	Rs. 50,000 (USD \$650) + shipping and taxes	~ Rs. 15,00,000	~ Rs. 1,30,000




WELSTAT



What is Potentiostat?


- Essential tool in electrochemical research
- It allows the user to apply a potential to a system (i.e. an electrochemical cell) and measure the resulting current, or vice versa
- It can do so while keeping the path where the current flows separate from the path where the potential is measured



Source: *Electrochemical Methods, Bard and Faulkner* Current axis output





Why do we need Potentiostat?



WELSTAT

- WELStat is a low cost, portable potentiostat
- Currently it can perform Cyclic voltammetry (techniques like DPV can be added)
- Used for detecting target DNA (amplified by our WELPCR)







Choosing Analog Front End IC (AFE)

We started with AD5941 to see whether it can be used for our work instead of EMSTAT Pico

AD5941



EMStat Pico




AD5941	EMStat Pico
₹ 700	₹ 45,000
Needs an external micro-controller	Micro-controller + Potentiostat
Not pre-programmed	Pre-programmed
Programming using embedded-C language	human readable scripting language for ease of programming

Parameter	Similarity
Excitation voltage range	-1.2 to +2 V
Max current	±3mA
ADC	16 bit
DAC	12bit
No. of Potentiostat	2
No. of channels	2


Although both EMSTAT Pico and AD5941 have almost similar specifications, the high-cost of EMSTAT Pico is mainly due to the well designed user-friendly software and accurately calibrated Potentiostat

Testing AD5941 for WELSTAT


(a) Testing AD5941 for our application using readymade Evaluation board



(b) Made PCB for AD5941 IC to check required minimal circuit for the IC




(c) Final PCB with AD5941 and micro-controller on same PCB



WELSTAT 2.0

- Both AFE AD5941 and Micro-Controller STM32 on same PCB
- Data can be read using USB or UART
- Status LED to indicate status of current running task



WELStat connected to mobile and mobile showing raw data


Challenges

1	Soldering QFN and QFP package ICs
2	Porting the Example program given by manufacture for STM32 controller
3	Electromagnetic interference (EMI)
4	Every time connector need to be soldered on electrode
5	Real time sensing while PCR reaction is in progress




Challenge: Electromagnetic Interference (EMI)





Electrode wire in contact with signal wire



Long wire captures signal of nearby digital devices



No wire




Challenges	Solution	Remark
Soldering QFN and QFP package ICs	Self-learnt soldering QFN and QFP package ICs	Successful
Porting the example program given by manufacturer for STM32 controller	Converted the program given for STM32F411 to STM32F103	Successful
Electromagnetic interference (EMI)	Care taken while PCB designing to minimize electrode track length as small as possible. No signal carrying track should be near the electrode tracks.	Successful
Every time connector needs to be soldered on electrode	USB male connector directly printed on electrode	Successful
Real-time electrochemical sensing of PCR amplicons	Need for further calibration of potentiostat and make it Noise-proof	Future scope



Electrochemical measurement using WELSTAT

Electrochemical measurements using cyclic voltammetry (CV) for different concentrations of methylene blue (MB)




Conclusion:
increase in MB
concentration results in an
increase in current




End point PCR using WELPCR+WELSTAT

- 503bp fragment of bacteriophage Phi-6 was amplified using WELPCR
- The amplicons were incubated with 50 μ M methylene blue for 1 hour at 4°C
- Samples were dispensed on PCB electrode and electrochemical measurements (CV) was performed using WELSTAT



Future Improvements

- Reduce the size of WELPCR device. The MOSFETs used in the current device require a big heat sink thereby occupying more space. These MOSFET, gate driver and heat sink can be replaced by VNH2SP30, which is an H-bridge motor driver IC to reduce the overall size of the device.



H-bridge using discrete MOSFET



VNH2SP30: H-bridge IC

- Reduce noise for WELSTAT to perform Real-time electrochemical PCR

Acknowledgement

- Wadhwani Electronics Lab, EE, IIT Bombay
- Prof. Siddharth Tallur, EE, IIT Bombay
- WEL Lab Staff Maheshwar Mangat, Mahesh Bhagnagare, Amit Shetye, Sadanand Sawant, V V Shaheen, Suraj Sarfare
- Mechanical and Tools Lab Staff Vishwas Kewale
- WEL RA Poulami Mandal, Bhagyashree Singh
- AIMS lab members Ruchira Nandeshwar, Avani Kulkarni, Shruti Ahuja, Rutuja Chalke



THANK YOU

WELPCR